

phototherapy **LAMP**S

006-PFL

**ROTARY PHOTOTHERAPY WITH
FLUORESCENT BULBS
ANGLE OF LAMP CAN BE ADJUSTED
HEIGHT OF LAMP MODULE CAN BE ADJUSTED
WHEELS CAN BE LOCKED**

***FOTOTERAPIA GIRATORIA CON
LAMPARAS FLUORESCENTES
EL ANGULO DE LA LAMPARA PUEDE SER AJUSTADO
LA ALTURA DE LA LAMPARA PUEDE SER AJUSTADA
LAS REUDAS PUEDEN SER BLOQUEADAS***



Phototherapy Lamps



STANDARD CONFIGURATION

Mobile model: The lamp (including lamp module, protective panel and timer), lock knob, and mobile stand (including height adjustment rod, lamp column and stand), and wheels.

TECHNICAL SPECIFICATIONS

Power Supply: AC110V-120V/60Hz

Maximum Power Output: 100W(XHZ-90)

Light source: blue fluorescent tube (XHZ-90)

Peak spectrum: 436nm(XHZ-90)

The highest total irradiance for bilirubin on the effective surface area:

$\geq 1.4 \text{mw/cm}^2$ (xhz-90 mobile)

Uniformity of the total irradiance for bilirubin:

> 0.5 (XHZ-90 mobile)

Wave length of blue light: 420-470nm

Lifetime of the light source: 2000hrs(blue fluorescent tube)

Adjust range for radiation angle of lamp module: 0-60° up from horizontal

Adjusting range for the height of lamp module:

1350mm -1650mm

Built-in timer range: 0-9999 hrs 59min

PACKAGE

Package: each unit is packed as one carton

Size: 95.5cm x 65.5cm

Gross Weight: 23Kg