

HOSPITAL BEDS

Camas para Hospital

B-500 Series

(SERIE B-500)



B-500E (Displayed)

Three Position Electric Bed
Camá Eléctrica de tres posiciones

4 Bed Panels
Camá de 4 Planos



B-500E TECHNICAL SPECIFICATIONS / ESPECIFICACIONES TECNICAS

Dimensions: Length 212cm, Width 100cm
Maximum/Minimum Height: 70cm/40cm
3 actuators, 1 control box, 1 handset
4 Bed Panels
Three Positions (Up-Down, Back, Knee)
Aluminum Alloy Bedside Rails
Six Pcs IV drip holes
Easily removable head and footboards
Central locking system (12.5cm castors)
Max weight capacity 175 kg (385 Lbs)
Overall net weight of bed 114 Kg (250 Lbs)
Voltage 110V / 60 Hz
2 urine hooks

<i>Dimensiones: Largo 212cm, Ancho 100cm</i>
<i>Altura Máxima/Mínima: 70cm/40cm</i>
<i>3 motores, 1 caja controladora, 1 Control de Mano</i>
<i>Cama de 4 Planos</i>
<i>Tres Posiciones (Arriba-Abajo, Espalda, Rodilla)</i>
<i>Barandas plegables de aluminio</i>
<i>Seis receptáculos para el atril porta suero</i>
<i>Cabecero superior e inferior (pies) desmontables</i>
<i>Sistema central de freno (Ruedas de 12.5cm de diámetro)</i>
<i>Capacidad máxima 175 Kg (385 Lbs)</i>
<i>Peso neto de la cama 114 Kg (250 Lbs)</i>
<i>Voltage 110V / 60 Hz</i>
<i>2 asas para bolsa de orina</i>

B-500M TECHNICAL SPECIFICATIONS / ESPECIFICACIONES TECNICAS

Dimensions: Length 212cm, Width 100cm
Maximum/Minimum Height: 70cm/40cm
4 Bed Panels
Three Positions (Up-Down, Back, Knee)
3 sets manual crank systems
Aluminum Alloy Bedside Rails
Six Pcs IV drip holes
Easily removable head and footboards
Central locking system (12.5cm castors)
Max weight capacity 175 kg (385 Lbs)
Overall net weight of bed 114 Kg (250 Lbs)
2 urine hooks

<i>Dimensiones: Largo 212cm, Ancho 100cm</i>
<i>Altura Máxima/Mínima: 70cm/40cm</i>
<i>Cama de 4 Planos</i>
<i>Tres Posiciones (Arriba-Abajo, Espalda, Rodilla)</i>
<i>3 Manivelas</i>
<i>Barandas plegables de aluminio</i>
<i>Seis receptáculos para el atril porta suero</i>
<i>Cabecero superior e inferior (pies) desmontables</i>
<i>Sistema central de freno (Ruedas de 12.5cm de diámetro)</i>
<i>Capacidad máxima 175 Kg (385 Lbs)</i>
<i>Peso neto de la cama 114 Kg (250 Lbs)</i>
<i>2 asas para bolsa de orina</i>

Contact us for additional information.
 Note: Specifications subject to change without notice.

Contáctenos para información adicional
 Nota: Especificaciones sujetas a cambios sin previo aviso.



6800 N.W. 77 Court, Miami, FL 33166
 Phone: 305-477-6331 • Fax: 305-477-5351
 sales@advanced-inst.com
www.advanced-inst.com