

UNIVERSAL OPERATING TABLE

MESA DE OPERACIONES UNIVERSAL

OT-30



Manual operating table
Mesa de operaciones Manual


SB-OT30
REV. 00
OCT. 2009

 **Advanced**

Manual Operating Table / Mesa de Operaciones Manual

TECHNICAL SPECIFICATIONS

ESPECIFICACIONES TÉCNICAS

Table Top/ <i>Tope de Mesa</i>	2000mmx500mm	
Height Without Pad/ <i>Altura sin colchon</i>	780mm-1020mm	
Trendelenburg-Reverse / <i>Trendelenburg-Invertido</i>	+22° / -22°	
Lateral Tilt (right)/ <i>Movimiento Lateral</i>	+18° / -18°	
Back Plate (up/down)/ <i>Espalda (arriba/Abajo)</i>	+75° / -15°	
Leg Plates (up/down)/ <i>Panel para las piernas</i>	+15° / -90°	
Leg Plates horizontal opening/ <i>Apertura para las piernas</i>	90°	
Flex/ <i>Quiebre de mesa</i>	225°	
