

ELECTRO HIDRAULIC SURGICAL TABLE

MESA DE CIRUGIA ELECTRO HIDRAULICA

OT-500

Standard Accessories:

Accesorios estandard

Anesthesia Supporter

Soporte para Anestesia

Hand Rest

Apoyo para la mano

Leg rest

Apoyo para la pierna

Chest & Waist Support

Soporte para el pecho y la cadera

Padding

Colchon

Body Fixed Belt

Cinturon fijo para el Cuerpo

Longitudinal Shift (optional)

Movimiento horizontal (opcional)



General Surgery



Kidney Surgery



Ortopedic Surgery



Urology



Gynecology

SB-OT500
REV. 00
OCT. 2009

Electro Hydraulic Surgical Table / Mesa de Cirugía Electro Hidraulica

TECHNICAL SPECIFICATIONS ESPECIFICACIONES TÉCNICAS

| | |
|--|-------------------------|
| Table Top/ <i>Tope de Mesa</i> | 2045mmx520mm |
| Height Without Pad/ <i>Altura sin colchon</i> | 710mm-1110mm |
| Trendelenburg-Reverse / <i>Trendelenburg-Invertido</i> | 22° / 18° |
| Head plate (up/down)/ <i>Soporte cabeza (arriba/abajo)</i> | 20° / 90° |
| Lateral Tilt (right)/ <i>Movimiento lateral</i> | 18° / 18° |
| Back Plate (up/down)/ <i>Espalda (arriba/abajo)</i> | +75° / -25° |
| Leg Plates (up/down)/ <i>Panel para las piernas (arriba/abajo)</i> | +20° / - 90° |
| Leg Plates horizontal extension/ <i>Apertura para las piernas</i> | 90° |
| Power/ <i>Voltage</i> | 110V, 60Hz / 220V, 50Hz |

Integrated the advanced ergonomic design and fine manufacturing arts, the OT500 can electrically adjust comfortably and easily to meet the requirements of varied operations. It has been approved by ISO 9001/13485 and CE, thus it becomes the best choice for present operating room.

