

# patient STRETCHER

## ST-1000 Plus



- *Base Durable Termoformada con Compartimiento para Accesorios*  
*Durable Thermoformed Base with Storage Tray*
- *Baranda Plegable*  
*Fold Down Side Rails*
- *Ruedas 8" para Fácil Transporte de Paciente*  
*8" Casters for Easy Patient Transport*

 **Advanced**

# ST-1000 PLUS

Patient Stretcher



## TECHNICAL SPECIFICATIONS

### STANDARD FEATURES

4" Mattress  
8" Casters  
27" Wide Patient Surface  
Fold Down Side Rails  
Side Rails 14 1/2" X 49"  
Four I.V. Poles Receptacles  
Stationary Foot Section  
Four Brake / Steer Pedals  
Dual Pedestal Hydraulics  
Dual Side Mounted Foot Pedal  
Control Hydraulics  
Hands Free Trendelenburg  
Reverse Trendelenburg  
One Detachable I.V. Pole

### PRODUCT SPECIFICATIONS

Overall Length: 83"  
Overall Width: 31"  
Mattress: 76" Length  
27" Width  
Height Adjustment Range:  
Min Height with Casters 23 1/2"  
Max Height with Casters 35" ( 89cm )  
Litter Position  
Trend. / Reverse Trend: +/- 18°  
Backrest Adjustment: + 90°

### CASTERS

8" With Central Locking One Tracking  
4 Brake / Steer Pedals  
Weight Capacity: 550 Lbs

## ESPECIFICACIONES TÉCNICAS

### CARACTERÍSTICAS ESTANDARD

Colchón 10.2 cm  
Ruedas 20.3 cm  
Superficie para Paciente 68.6 cm  
Barandas Plegables  
Barandas 36.9 cm X 124.5 cm  
4 Receptáculos para paral de Suero  
Sección de los Pies Fija  
Cuatro Pedales de Freno  
Pedestal Hidráulico Doble  
Pedales en Ambos Lados  
Control Hidráulico  
Trendelemburg Manos Libre  
Trendelemburg Invertido  
Un Paral de Suero Desmontable

### ESPECIFICACIONES

Largo Completo 211 cm  
Ancho Completo 79 cm  
Largo del Colchón 193 cm  
Largo del Colchón 68.6 cm  
Altura Ajustable  
Altura Mínima con Ruedas 59.7 cm  
Altura Máxima con Ruedas 89 cm  
Posiciones  
Trendelemburg / Trendelemburg Invertido: +/- 18°  
Ajuste de la Espalda + 90°

### RUEDAS

Ruedas de 8" con Sistema Central de Freno  
4 Ruedas con Freno / Pedal de Dirección  
250 Kg Capacidad Máxima

Contact us for additional information.  
Contáctenos para información adicional.

NOTE: Specifications subject to change without prior notice.  
NOTA: Especificaciones sujetas a cambio sin previo aviso.

6800 N.W. 77 Court • Miami, Florida 33166 • USA  
Phone: 305-477-6331 • Fax: 305-477-5351  
info@advanced-inst.com  
[www.advanced-inst.com](http://www.advanced-inst.com)