


DxWorks

DR WORKSTATION



DxWorks

Direct Digital with Improved Workflow

NMI

Enhanced Workflow

The DxWorks console provides workflow improvement by automating all aspects of radiographic image acquisition, without eliminating manual operator control. Patient demographics are easily generated or transferred from DICOM worklist. Exposure and image processing parameters are pre-selected and modified by the operator.

Full resolution and processed images can be viewed and manipulated on the console prior to be sent to PACS for diagnosis and archive. DxWorks offers a complete set of image tools for full, professional image control.

- DR Image Acquisition
- Automated Detector Calibration
- Automated Image Processing and Fast Review
- Easy to Learn & Use
- Full DICOM Connectivity
- Optional Direct Generator Control
- Organ Specific Image Enhancement (OSIE™)
- Built-in Stitching Function
- True 14-bit Image Processing
- Image Disk Archive (DICOM DIR)



Backed by extensive experience in medical diagnostic equipment and the optical electronics industry, NMI is able to leverage its global research and development footprint to bridge the healthcare and IT industries. NMI is positioned to push new fronts in both technology and clinical applications, providing enhanced diagnostic capabilities to imaging centers through innovative and flexible products.